

****VaDia readme - Manager ****

This file "readmeVaDia - Manager " is the Vadia Manager.exe readme file.

CHANGE LOGS:

ver. 2.01 – 01.01.2014 file: VaDia Manager.exe
changes:

- Improvement in battery level indicator.
- Removed menus: Setup, View.

ver. 1.15 – 25.02.2013 file: VaDia Manager.exe
changes:

- added troubleshooting messages.
- button "Status" refreshes the USB connection
- more safety user interface

ver. 1.14 – file: VaDia Manager.exe
changes:

- not released.

ver. 1.13 - 15.12.2011 file: VaDia Manager.exe
changes:

- Automatic prefix 'VD'.

ver. 1.12 - 13.11.2011 file: VaDia Manager.exe
changes:

- Faster data transmission.
- New modes of VaDia. STOP, ERROR, UNDEFINED

ver. 1.11 - 17.10.2011 file: VaDia Manager.exe
changes:

- New text for warning box before deleting data and recalibrating Vadia.

ver. 1.10 - 03.09.2011 file: VaDia Manager.exe
changes:

- Warning box before deleting data and recalibrating Vadia.

ver. 1.09 - 01.01.2011 file: VaDia Manager.exe
changes:

- initial version